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1. PRODUCT DESCRIPTION

The Fergus-brand suture line consists of the CordFit, CordFit A, TapeFit, TapeFit A and TapeFit Max models. They are surgical, sterile, non-absorbable sutures manufactured with braided multifilaments made of ultra-high-molecular-weight polyethylene (UHMWPE) or a mixture of UHMWPE and polypropylene or polyamide (nylon).

The sutures are designed for applications that involve approximating or joining soft tissue.

CordFit and CordFit A suture models are thread-type sutures, and TapeFit and TapeFit A sutures are flat or tape-type sutures. TapeFit Max is a combination of both previous models: a thread-type suture on the ends and a tape-type suture in the center.

Fergus-brand sutures are supplied in preset lengths and with stiffened ends for easier handling and threading.

This product line features the following sutures:

- UHMWPE sutures: in white, white-blue, and blue. The Cobalt Chromite Blue-Green Spinel pigment is used for the blue threads. White threads are not pigmented.
- UHMWPE and polypropylene sutures: partially pigmented following a specific pattern. The pigment used for the blue part of the white-blue suture is phthalocyanine. White threads are not pigmented.
- UHMWPE and polyamide sutures: partially pigmented following a specific pattern. The pigment used for the black part of the white-black suture is hematein. White threads are not pigmented.

CordFit sutures come in nine models, with different sizes and colors, as follows:

- White: #0, #2 and #2-0 UHMWPE sutures, 38" (96 cm) in length.
- White/blue: #0, #2 and #2-0 UHMWPE sutures, 38" (96 cm) in length.
- Blue: #0, #2 and #2-0 UHMWPE sutures, 38" (96 cm) in length.

CordFit A sutures come in four models, with different sizes and colors, as follows:

- White/blue: #2 and #2-0 UHMWPE and polypropylene sutures, 38" (96 cm) in length.
- White/black: #2 and #2-0 UHMWPE and polyamide sutures, 38" (96 cm) in length.

TapeFit sutures come in nine models, with different sizes and colors, as follows:

- White: 1.5 mm, 2 mm and 3 mm UHMWPE sutures, 38" (96 cm) in length.
- White/blue: 1.5 mm, 2 mm and 3 mm UHMWPE sutures, 38" (96 cm) in length.
- Blue: 1.5 mm, 2 mm and 3 mm UHMWPE sutures, 38" (96 cm) in length.

TapeFit A sutures come in four models, with different sizes and colors, as follows:

- White/blue: 1.5 mm, 2 mm and 3 mm UHMWPE and polypropylene sutures, 38" (96 cm) in length.
- White/black: 1.5 mm UHMWPE and polyamide sutures, 38" (96 cm) in length.

Finally, TapeFit Max sutures come in two models, with different colors, as follows:

- White/blue: 2 mm UHMWPE and polypropylene sutures, 254-254-254 mm (10"-10"-10") in length.
- White/black: 2 mm UHMWPE and polyamide sutures, 254-254-254 mm (10"-10"-10") in length.

2. PRODUCT PRESENTATION

The CordFit, CordFit A, TapeFit, TapeFit A and TapeFit Max Fergus sutures are supplied in 6- or 12-unit boxes. Each product unit is packaged individually. Check the individual suture labels for the model, material, color, size and length. All components are provided sterile and ready for use.

3. INDICATIONS FOR USE

Fergus-brand sutures are intended for applications that involve approximating or joining soft tissue, and can be used as components in procedures that include repairs using allograft or autograft tissue in arthroscopic surgery, minimally-invasive orthopedic surgery and/or open orthopedic surgery.

4. CONTRAINDICATIONS

The product should not be implanted in the following cases:

- Quality or pathologies of the tissue to be reconstructed which may affect the safety of the attachment by means of suture.
- Limited blood irrigation and/or previous infection which may delay healing and/or consolidation.
- Heart surgery, wound closing or surgical procedures not mentioned in the "Indications for Use" section.

5. WARNINGS

- Read the instructions for use carefully and completely before use.
- The product must be used only by a medical professional trained in the surgical techniques indicated for placing these implants.
- The law allows this product to be sold only to physicians or people authorized by a physician.
- DO NOT use the product if the packaging is open or damaged, since product sterility may be lost.
- The product has been designed for SINGLE use only. Therefore, IT MUST NOT BE REUSED, REPROCESSED OR RSTERILIZED, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative or governmental requirements.
- Do not use the product after the expiration date that appears on the label. Product performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- Do not expose the suture to a direct heat source.
- Make sure that the product is disposed of as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- The product must be stored in the original packaging.
- Patient sensitivity to product materials must be considered before implantation.
- All surgical instruments used in the surgical procedure to install this implant must be sterile and suitable for use in aseptic environments.
- The fixation provided by the suture must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must

be followed carefully in order to prevent overstressing the implant.

- Patients will receive detailed instructions about the use and limitations of the implant.
- If the product is to be explanted, the risks of the possibility of a second procedure on the patient must be considered. Product extraction must be accompanied by an appropriate post-operative approach.
- Sutures that have been opened or introduced into the patient must be disposed of, even if they have not been used or handled, since this could increase the risk of inappropriate sterilization and cross-contamination.
- Bio-hazardous waste, such as explanted implants, sutures, needles and contaminated surgical materials, must be disposed of safely and as per site policy.
- Sutures are non-absorbable and may become encapsulated in the adjacent connective tissue.
- The suture has not been found or proven to be safe and effective for use as an artificial tendon or ligament.
- The manufacturer shall not be held liable for faults or damages caused when these implants are combined with products not recommended by the manufacturer.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, which may affect patient safety and implant performance.

6. PRECAUTIONS

- Surgeons must follow their professional judgment to determine the appropriate implant size and model based on the specific indication and the history of the patient.
- Surgeons are advised to review the product-specific surgical technique before the procedure.
- Before using the product, examine it for damage. Do not use if damaged.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper product selection and installation, are essential for satisfactory use of the implant.
- Sutures may break if handled or controlled using sharp or cutting instruments.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.

7. SUPPLY AND STORAGE

- The sutures are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- In order to maintain product quality and effectiveness, they are supplied in a cardboard container and a double Tyvek pouch. After the sterilization, the product is protected during transportation and storage by placing it inside a cardboard box (external packaging). Labels may be added for identification.
- Keep the product away from any sources of heat or humidity.
- Room temperature storage is advised.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- Implant shelf life is 3 years.

8. MAGNETIC RESONANCE (MR) ENVIRONMENTS

- MR Safety and Compatibility

Sutures are made of ultra-high-molecular-weight polyethylene (UHMWPE) or a mixture of UHMWPE with polypropylene or polyamide (nylon), which makes them MR safe.

9. POSSIBLE COMPLICATIONS

All surgical procedures involve a certain risk. Some of the complications that could arise are:

- Deep and superficial infections.
- Reactions to foreign bodies.
- Pain on the site of the incision or the procedure.
- Wound dehiscence.
- Edema and erythema on the site of the wound.
- Gradual suture encapsulation in fibrous connective tissue.
- Tissue damage.
- Incomplete or inappropriate lesion correction, or lesion recurrence.
- Possible allergies and other types of reactions to implant materials.

10. PATIENT INFORMATION

Surgeons are responsible for informing patients and their representatives before surgery about the possible complications related to the implantation (see the "Possible Complications" section).

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated with the implantation and use of the product.

Post-operative care is important. Patients will receive instructions about the limitations required by the implant and will be warned about the necessary precautions they must apply for healing and lesion resolution.

Patients should be warned that if they decide to explant the product, they should consider the risks of undergoing a second surgical procedure. Product removal should be accompanied by an appropriate surgical and post-surgical approach.

Performing additional surgeries does not guarantee a complete repair of the complications.

11. SURGICAL TECHNIQUE

Surgeons must be trained in the use accepted techniques for surgical knots or knotless repairs, suture numbers and configurations for the surgical procedures mentioned in the "Indications for Use" section, the types of lesions to be repaired, the clinical history of the patient and other criteria that the professional considers important to achieve a satisfactory use of the product during the surgical procedure.

This product is often used in combination with other medical products following the surgical technique associated to those products.

After use, discard the product and its packaging per local hospital, administrative or governmental procedures.

12. MATERIAL SPECIFICATIONS

CordFit y TapeFit suture models are made with ultra-high-molecular-weight polyethylene (UHMWPE). The pigment used is Cobalt Chromite Blue-Green Spinel for the blue color. White fibers in the suture are not pigmented.

CordFit A, TapeFit A and TapeFit Max suture models are manufactured with a mixture of UHMWPE and a smaller percentage of polypropylene for white-blue suture, or a mixture of UHMWPE and a smaller percentage of polyamide (nylon) for white-black suture. The blue pigment used for these sutures is phthalocyanine, while the black pigment used is hematein. White fibers in the suture are not pigmented.

Check the packaging labels for information about the materials.

13. LEARN MORE

To learn more about this product, visit our website:
www.fergusmedical.com

14. SYMBOLS USED IN LABELS

| | |
|---|--|
|  | CATALOG NUMBER |
|  | BATCH NUMBER |
|  | SERIAL NUMBER |
|  | PRECAUTION |
|  | EXPIRATION DATE |
|  | DO NOT REUSE |
|  | MANUFACTURING DATE |
|  | STERILE. STERILIZATION METHOD: ETHYLENE OXIDE |
|  | DOUBLE STERILE BARRIER SYSTEM |
|  | MANUFACTURER |
|  | CHECK INSTRUCTIONS FOR USE |
|  | DO NOT USE IF PACKAGING IS DAMAGED |
|  | DO NOT RE-STERILIZE |
|  | KEEP AWAY FROM SUNLIGHT |
|  | KEEP DRY |
|  | MR-SAFE |

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