# **Knotix PK**

# **QuickFit PK**

# **Knotix SPD PK**

# **SpinFit PK**

---

Anclas de PEEK

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Instructions for use  
ENGLISH

M-351-00148 (Rev.07) / 25.03.2025

# Fergus

Fergus

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**MANTER SECO**

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## 1. DEVICE DESCRIPTION

The Fergus-brand PEEK anchors include the Knotix PK, Knotix SPD PL, QuickFit PK and SpinFit PK models. These are PEEK devices that are designed to be threaded or impacted and feature eyelets to allow for the passage of suture. They are supplied sterile and inside a disposable insertion device to aid in the placement of the implant.

Some PEEK anchors include non-absorbable sutures, while others do not, and feature a suture threader to help threading a suture afterwards.

The PEEK anchors are designed to attach soft tissue to the bone.

The **Knotix PK** model features a threaded anchor and an internal eyelet for two CordFit #2 sutures (USP standard). Sutures may come in blue, white or white-blue colors. The four tips of these two sutures are connected to the suture retention device located in the handle of the insertion device. This model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 5.5 mm diameter, 14.7 mm length.

Figure 1: Knotix PK PEEK anchor.

The **Knotix SPD PK** model comes in two variants:

- an anchor with spikes for impact implantation and a CordFit #2 (USP standard) non-absorbable suture threaded in the external eyelet of the anchor. Sutures may come in blue, white or white-blue colors. The 2 tips of the suture are connected to the suture retention device located in the handle of the insertion device. The Knotix SPD PK model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 3.0 mm diameter, 12.4 mm length.
- an anchor with spikes for impact implantation and an internal eyelet for two non-absorbable CordFit #2 sutures (USP standard). Sutures may come in blue, white or white-blue colors. The four tips of these two sutures are connected to the suture retention device located in the handle of the insertion device. The Knotix SPD PK model requires that the suture is to be connected by tying knots. This anchor is available in the following dimensions: 3.0 mm diameter, 12.2 mm length.

Figure 2: Knotix SPD PK PEEK anchor  
(Top: one suture. Bottom: two sutures)

The **QuickFit PK** PEEK anchor is a two-part implant: an anchor with spikes for impact implantation and the eyelet to pass the non-absorbable suture. The suture is not included and must be threaded by the physician. This model includes a suture threader built into the eyelet to aid in threading the suture. The eyelet is attached to the insertion device by means of a small thread in the rod. During implantation, the anchor is connected to the eyelet, which allows for the suture to remain in place due to an interference phenomenon between the bone tissue and the implant, with no knots required. This model is available in 4.5 mm diameter and 24 mm length or 2.9 mm diameter and 15.5 mm length. Dimensions are identified in product labels.

**NOTE:** The length stated equals the sum of the lengths of the anchor and the eyelet.

Figure 3: QuickFit PK PEEK anchor

The **SpinFit PK** PEEK anchor is a two-part implant: a threaded anchor and an eyelet to pass the suture. The SpinFit PK anchor is available with or without non-absorbable sutures in the following dimensions: 4.7 mm diameter and 19.1 mm length, or 3.5 mm diameter and 15.8 mm length. Dimensions are identified in product labels.

**NOTE:** The length stated equals the sum of the lengths of the anchor and the eyelet.

The anchor eyelet is distal, separated from the anchor and attached to the insertion device by means of an eyelet retention suture, for the 4.7-mm-diameter anchor, or a thread, for the 3.5-mm-diameter anchor. The SpinFit PK anchor model that does not include a suture, comes with a suture threader to aid the physician in passing the sutures through the eyelet. If the model includes a suture, it is a TapeFit Max model suture, which is a combination of a thread-type suture on the ends (#5, USP standard) and a 2-mm tape-type suture in the center. It has a 30" length and it is supplied in white-blue or white-black color. During implantation, the anchor moves and is attached to the eyelet when they come into contact. This causes the suture to become fixed in place by means of an interference phenomenon between the bone tissue and the implant thread, with no knots required.

Figure 4: SpinFit PK PEEK anchor with suture

Figure 5: SpinFit PK PEEK anchor without suture

## 2. PRODUCT PRESENTATION

Contents:

- **Knotix PK models:** one threaded PEEK anchor with two CordFit #2 implantable sutures and one insertion device.
- **Knotix SPD PK models:**
  - One spiked PEEK anchor with one CordFit #2 implantable suture and one insertion device.
  - One spiked PEEK anchor with two CordFit #2 implantable sutures and one insertion device.
- **QuickFit PK models:** one spiked PEEK anchor, one insertion device and one suture threader.
- **SpinFit PK model:**
  - One threaded PEEK anchor, one insertion device and one suture threader.
  - One threaded PEEK anchor with one TapeFit Max suture and one insertion device.

Check the individual labels of the device to find the model, anchor size, materials, whether the suture or the insertion device are included or not, and the color of the suture. All components are provided sterile and ready for use.

The availability of models and indications for use may vary from one market to another, which is subject to regulatory approvals in each country.

Please, check the Indications for use section to find the indications authorized for each country.

If you have any questions, please direct them to the commercial representative in your country.

## 3. INDICATIONS FOR USE

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK PEEK anchors are intended for attaching soft tissue to the bone. They are designed to be used in the shoulder, hip, knee, and ankle.

**NOTE:** In Mexico, they are authorized to be used in the shoulder and ankle only.

## 4. CONTRAINDICATIONS

The device should not be implanted in the following cases:

- Insufficient bone quality or quantity for attachment. The physician must evaluate carefully the bone quality of patients with bone immaturity before starting with the orthopedic surgical procedure, and must avoid deviating, interrupting or altering the growth cartilage.
- Quality or pathologic conditions of the soft tissue to be fixed which could affect the safety of fixation with the use of sutures.
- Limited blood irrigation or previous infection which may delay healing.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions of their physical activity.

- Surgical procedures not mentioned in the “Indications for Use” section.
- Inappropriate neuromuscular conditions.

## 5. WARNINGS

- Read the instructions for use carefully and completely before use.
- The device must be used only by a medical professional trained in the surgical techniques indicated for placing these implants.
- The law allows this device to be sold only to physicians or people authorized by a physician.
- DO NOT use the product if the packaging is open or damaged, since product sterility may be lost.
- The products have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative and/or governmental requirements.
- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- Make sure that the product is disposed of as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- The product must be stored in the original packaging.
- Patient sensitivity to device materials must be considered before implantation.
- All surgical instruments used in the procedure to install this implant must be sterile and suitable for use in aseptic environments.
- Product presentation must not be altered in a way different from what is stated in the surgical technique, since this could affect implant performance.
- If the implant detaches from the insertion device, DO NOT REATTACH IT, since this could impair the effectiveness and increase the risk of sterility loss.
- The fixation provided by the implant must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing the product and the bone.
- Patients will receive detailed instructions about the use and limitations of the product.
- If the decision is made to explant the product, the risks of undergoing a second surgical intervention should be considered. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.
- Opened devices and implants must be disposed of, even if they have not been used or handled, since this could increase the risk of inappropriate sterilization and cross-contamination.

- Bio-hazardous waste, such as explants, sutures, instruments and contaminated surgical materials, must be disposed of safely and as per site policy.
- The anchor and/or the insertion device may break when the insertion site is not properly prepared with the appropriate instrument before using the implant.
- Do not use the product if any faults are detected.
- Do not use any products or sutures that are not recommended by the manufacturer, since implant performance for the repair cannot be guaranteed in this case.
- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the manufacturer.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, which may affect patient safety and implant performance.

## 6. PRECAUTIONS

- Surgeons must follow their professional judgment to determine the appropriate implant model, quantity and size based on the specific indication and the history of the patient.
- The physiological measurements of the patient limit the implant and insertion device sizes that may be used. The professional must select the type and size that adapts better to the needs of the patient.
- Surgeons are advised to review the product-specific surgical technique before the procedure.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper device selection and installation, are essential for satisfactory use of the implant.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.
- PEEK anchor implantation requires preparation in the insertion site.
- The surgeon must consider the manufacturer-recommended installation systems in order to implant the product properly.
- The surgeon must consider the anchor size to select the punch, drill bit or bone drill, since an incomplete insertion of the anchor may cause a deficient performance of the anchor or an inappropriate implantation.
- If the bone drill or punch is not inserted to the appropriate depth, the anchor may become damaged or broken.
- In order to prevent damage on the anchor and/or insertion device, the surgeon must insert the product in the same orientation as the perforation in the bone.
- Using excessive force during insertion may cause faults in the suture anchor or the insertion device.

- Sutures may break if handled or controlled using sharp or cutting instruments.
- The decision of removing the implant must be made considering the potential risk involved in a second surgical procedure for the patient. After the explantation, appropriate post-operative controls must be carried out.

## 7. SUPPLY AND STORAGE

- PEEK Anchors are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- In order to maintain product quality and effectiveness, it is marketed in a double Tyvek pouch (primary packaging). After sterilization, the product is protected during transport and storage by placing it inside a cardboard box (external packaging). Additionally, the product is supplied with a plastic protective wrap covering the implant, which is placed to prevent Tyvek pouch damage. The product is identified with labels.
- Keep the product away from any sources of heat or humidity.
- Room temperature storage is advised.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- Product shelf life is 3 years.

## 8. MAGNETIC RESONANCE (MR) ENVIRONMENTS

- MR Safety and Compatibility

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK anchors are manufactured with polyether ether ketone (PEEK), which makes them MR safe.

Sutures are made with UHMWPE or a mixture of UHMWPE with polypropylene or polyamide (nylon), which makes them MR safe.

## 9. POSSIBLE COMPLICATIONS

All surgical procedures involve a certain risk. Some of the complications that could arise are:

- Deep and superficial infections
- Atypical reactions to foreign bodies
- Pain on the site of the incision or the procedure
- Implant loosening and/or detachment
- Bone lesions
- Incomplete or inappropriate lesion correction, or lesion recurrence.

## 10. PATIENT INFORMATION

The surgeon is responsible for informing patients and their representatives before surgery about the possible complications related to the implantation of the device (see the "Possible Complications" section).

The manufacturer and its distributors delegate to the surgeon the responsibility of informing patients about the advantages and possible risks associated with the implantation and use of the product.

Post-operative care is important. Patients will receive instructions from the physician about the limitations required by the implant, and will be warned about the necessary precautions they must apply for healing.

Patients should be warned that if they decide to explant the product, they should consider the risks of undergoing a second surgical procedure. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.

Performing additional surgeries does not guarantee a complete repair of the complications.

## 11. SURGICAL TECHNIQUE

The surgeon must select the appropriate approach based on the repair to be performed. The surgeon may opt for an open, minimally-invasive, or arthroscopic approach, as preferred.

Additionally, the surgeon must select the appropriate number of implants and fixation configuration based on the type of lesion, the anatomical region, the medical history of the patient and other criteria that the professional may deem important for a satisfactory use of the product during the surgical procedure.

The following surgical technique instructions refer only to device implantation, not considering other procedures that are critical to perform the surgical procedure, such as patient preparation, working portals, procedure finalization, etc.

Some interventions associated to each PEEK anchor model are suggested below. However, other procedures may be performed if the physician deems it appropriate as long as the indications for use are followed:

### ***Knotix PK***

Knotix PK is recommended for:

Shoulder	<ul style="list-style-type: none"><li>• Rotator cuff repair</li><li>• Proximal biceps tenodesis</li></ul>
Knee	<ul style="list-style-type: none"><li>• Lateral extra-articular tenodesis (LET)</li><li>• Posterolateral complex reconstruction</li><li>• Posteromedial complex reconstruction</li><li>• Medial patellofemoral ligament (MPFL) reconstruction</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Lateral instability repair</li><li>• Medial instability repair</li><li>• Achilles tendon repair</li></ul>

## Steps for proper Knotix PK anchor implantation:

**1) Shoulder:** Identify the anchor insertion site and place the punch for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN001047). Impact the end of the handle of the punch successively until the first laser marking of the instrument has been inserted completely in order to create the hole in the bone.

**ATTENTION:** for high bone mineral density, we advise using the drill bit for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN081047) and, afterwards, the threaded punch for 5.5-mm Knotix (FG-IN001155) instead of the punch for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN001047).

**Knee and ankle:** Identify the insertion site for the anchor, attach the drill bit for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN081047) to the drill, and perforate until reaching the reference point. Afterwards, using the threaded punch for 5.5-mm Knotix (FG-IN001155) create the imprint to facilitate the installation of the anchor until the entire thread of the instrument has been inserted.

**2)** Position the anchor in the implantation site.

**3)** Turn the handle of the insertion device clockwise until the transversal laser marking located on the distal end of the capillary (of the insertion device) has been completely inserted. Applying a slight forward (proximal-distal) pressure is advised while the handle is turned in order to aid in the implantation.

**ATTENTION:** the longitudinal capillary laser marking indicates the direction towards which the suture comes out. Assess the orientation of the laser marking based on the desired repair configuration.

**4)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor.

**5)** For arthroscopic approach: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached; tie the knot, and then cut the excess suture (if appropriate).

For open approach: the suture must be threaded in a needle in order to be passed through the tissue, and attached to the anchoring point. Finish the attachment by tying surgical knots.

**ATTENTION:** the needle is not included in the product.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

## **Knotix SPD PK**

The Knotix SPD PK anchor is recommended for:

Shoulder	<ul style="list-style-type: none"><li>• Glenoid labral repair</li><li>• Hill-Sachs lesion repair</li></ul>
Hip	<ul style="list-style-type: none"><li>• Acetabular labral repair</li></ul>
Knee	<ul style="list-style-type: none"><li>• Medial patellofemoral ligament (MPFL) reconstruction</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Lateral instability repair</li><li>• Medial instability repair</li></ul>

Steps for proper Knotix SPD PK anchor implantation (with one or two sutures):

- 1)** Identify the implantation site for the anchor and position the 3.2-mm implantation guide (FG-IN013000) in the selected insertion site.
- 2)** Insert the drill bit for 3.0-mm Knotix SPD (FG-IN082030) through the guide and drill in the insertion point until the drill bit reaches its stop on the guide. Then, remove the drill bit while keeping the guide in place.
- 3)** Insert the anchor through the guide and impact the handle successively until the transversal laser marking, located on the distal end of the capillary of the insertion device, reaches the edge of the implantation guide.  
**ATTENTION:** the longitudinal capillary laser marking indicates the direction towards which the suture comes out. Assess the orientation of the laser marking based on the desired repair configuration.
- 4)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor and the implantation guide.
- 5)** For arthroscopic approach: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached; tie the knot, and then cut the excess suture.

For open approach: the suture must be threaded in a needle in order to be passed through the tissue, and attached to the anchoring point. Finish the attachment by tying surgical knots.

**ATTENTION:** the needle is not included in the product.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

## SpinFit PK (4.7 mm in diameter)

The 4.7-mm-diameter SpinFit PK anchor is recommended for:

Shoulder	<ul style="list-style-type: none"><li>• Rotator cuff repair</li><li>• Proximal biceps tenodesis</li></ul>
Knee	<ul style="list-style-type: none"><li>• Lateral extra-articular tenodesis (LET)</li><li>• Posterolateral complex reconstruction</li><li>• Posteromedial complex reconstruction</li><li>• Medial patellofemoral ligament (MPFL) reconstruction</li><li>• Meniscal root repair</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Augmentation of instability repair</li><li>• Achilles tendon repair</li></ul>

Steps for proper 4.7-mm SpinFit PK anchor implantation:

**1) Shoulder:** Identify the anchor insertion site and place the punch for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN001047). Impact the end of the handle of the punch successively until the first laser line of the instrument has been inserted completely in order to create the hole in the bone.


**ATTENTION:** for high bone mineral density, we advise using the drill bit for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN081047) and, afterwards, the threaded punch for 4.7-mm SpinFit (FG-IN001147) instead of the punch for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN001047).

**Knee and ankle:** Identify the insertion site for the anchor, attach the drill bit for 5.5-mm Knotix and 4.7-mm SpinFit (FG-IN081047) to the drill, and perforate until reaching the reference point. Then, using the threaded punch for 4.7-mm SpinFit (FG-IN001147) create the imprint to facilitate the installation of the anchor until the entire thread of the instrument has been inserted.

**2)** Remove the 4.75-mm SpinFit PK punch and position the SpinFit PK anchor inside the hole created in the bone until the laser marking in the anchor capillary is not arthroscopically visible.

For SpinFit PK without TapeFit Max suture: before this step, the suture holding the tissue must be passed through the anchor eyelet using the disposable threader connected to the capillary of the insertion device.

**ATTENTION:** Check for the desired suture tension, since, from the insertion stage on, the suture will be secured, and you will be unable to readjust the attachment tension. To achieve this, leave a free suture length equaling the length of the anchor, i.e. the distance between the tip of the capillary and the laser marking.

Check that the distal part of the anchor comes into contact with the bone before introducing the body of the implant in the prepared bone perforation.

Uncouple the eyelet retention suture from the proximal end of the handle.

**3)** Turn the proximal end of the handle clockwise holding the insertion device to be able to move forward towards the inside of the hole in the bone. Turn until the implantable anchor is not visible and is inserted completely in the hole in the bone. Remove the eyelet retention sutures.

**4)** Release the TapeFit Max suture from the handle tabs and remove the insertion device from the anchor. For SpinFit PK without TapeFit Max suture: cut the excess suture.

**5)** For repairs: use the suture passing instrument to pass the sutures of the anchor through the tissue to be attached.

**6)** Cut the excess suture.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures

### **SpinFit PK** (3.5 mm in diameter)

The 3.5-mm-diameter SpinFit PK anchor is recommended for:

Shoulder	<ul style="list-style-type: none"><li>• Glenoid labral repair</li><li>• Hill-Sachs lesion repair</li></ul>
Knee	<ul style="list-style-type: none"><li>• Lateral extra-articular tenodesis (LET)</li><li>• Posterolateral complex reconstruction</li><li>• Posteromedial complex reconstruction</li><li>• Medial patellofemoral ligament (MPFL) reconstruction</li><li>• Meniscal root repair</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Augmentation of instability repair</li><li>• Achilles tendon repair</li></ul>

Steps for proper 3.5-mm SpinFit PK anchor implantation:

**1)** Identify the insertion site for the anchor and position a guide with an inner diameter of, at least, 3.2 mm, in the insertion site.

**2)** Insert a 2.9 mm drill bit through the inside of the guide and drill a hole 15.8 mm in depth in the insertion point. Remove the guide and the drill bit.

**3)** Thread the suture holding the tissue through the eyelet of the SpinFit PK anchor and hold the sutures in the suture retention area of the handle.

**ATTENTION:** Check for the desired suture tension, since, from the insertion stage on, the suture will be secured, and you will be unable to readjust the attachment tension. To achieve this, leave a free suture length equaling the length of the anchor, i.e. the distance between the tip of the capillary and the laser marking.

**4)** Release the sutures from the suture retention area of the handle and remove the insertion device from the anchor by rotating it counter-clockwise until it detaches from the eyelet. Cut the excess suture.

**5)** Turn the proximal end of the handle clockwise holding the insertion device to be able to move forward towards the inside of the hole in the bone. Turn until the implantable anchor is not visible and is inserted completely in the hole in the bone.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **QuickFit PK** (2.9 mm in diameter)

The 2.9-mm-diameter QuickFit PK anchor is intended for:

Shoulder	<ul style="list-style-type: none"><li>• Glenoid labral repair</li><li>• Hill-Sachs lesion repair</li></ul>
Hip	<ul style="list-style-type: none"><li>• Acetabular labral repair</li></ul>
Knee	<ul style="list-style-type: none"><li>• Medial patellofemoral ligament (MPFL) reconstruction</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Augmentation of instability repair</li></ul>


Steps for proper 2.9-mm QuickFit PK anchor implantation:

**1)** Identify the anchor insertion site and position the 3.2-mm implantation guide (FG-IN013000) in the insertion site.

**2)** Insert a drill bit for 2.9-mm QuickFit (FG-IN082029) through the inside of the guide and drill a hole in the insertion point. Remove the guide and the drill bit.

**3)** Thread the suture holding the tissue through the eyelet of the QuickFit PK anchor.

**4)** Adjust the suture tension, since, after the placement, the suture will be secured, and you will be unable to readjust the attachment tension. To achieve this, leave a free suture length equaling the length of the anchor, i.e. the distance between the tip of the capillary and the first laser marking.

**5)** Position the anchor over the perforation made and impact the handle successively to move the end of the inner rod (where the eyelet is located) forward into the perforation. Impacts should continue until the first laser marking (located in the rod) has been inserted completely in the perforation.

**6)** Remove and discard the safety clip located in the proximal end of the handle.

**7)** Impact the handle successively to advance it forward and CE the final implantation of the anchor. Impacts should continue until the second laser marking (located in the external capillary) has been inserted completely in the perforation.

**ATTENTION:** In order to prevent possible damage in the surrounding tissue or implant detachment, it is important that you ensure a complete insertion of the anchor. To achieve this, the impacts must be performed ensuring that the laser markings are not visible after each of the steps above.

**8)** Remove the insertion device from the anchor by rotating it counterclockwise until it detaches from the eyelet. Cut the excess suture.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

### **QuickFit PK** (4.5 mm in diameter)

The 4.5-mm-diameter QuickFit PK anchor is recommended for:


Shoulder	<ul style="list-style-type: none"><li>• Rotator cuff repair</li></ul>
Ankle	<ul style="list-style-type: none"><li>• Augmentation of instability repair</li><li>• Achilles tendon repair</li></ul>

Steps for proper 4.5-mm QuickFit PK anchor implantation:

**1)** Identify the anchor insertion site and position the 4.5-mm QuickFit PK punch. Impact the end of the handle of the 4.5-mm QuickFit PK punch successively until the second laser-marked line of the capillary of the instrument has been inserted completely in order to create the hole in the bone.

**2)** Thread the suture holding the tissue through the eyelet of the QuickFit PK anchor.

**3)** Adjust the suture tension, since, after the placement, the suture will be secured, and you will be unable to readjust the attachment tension. To achieve this, leave a free suture length equaling the length of the anchor, i.e. the distance between the tip of the capillary and the laser marking.

**4)** Position the anchor over the perforation made and impact the handle successively to move the end of the inner rod (where the eyelet is located) forward into the perforation. Impacts should continue until the first laser marking (located in the rod) has been inserted completely in the perforation.

**5)** Remove and discard the safety clip located in the proximal end of the handle.

**6)** Impact the handle successively to advance it forward and CE the final implantation of the anchor. Impacts should continue until the second laser marking (located in the external capillary) has been inserted completely in the perforation.

**ATTENTION:** In order to prevent possible damage in the surrounding tissue or implant detachment, it is important that you ensure a complete insertion of the anchor. To achieve this, the impacts must be performed ensuring that the laser markings are not visible after each of the steps above.

**7)** Remove the insertion device from the anchor by rotating it counterclockwise until it detaches from the eyelet. Cut the excess suture.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

## 12. MATERIAL SPECIFICATIONS

The Knotix PK, Knotix SPD PK, QuickFit PK and SpinFit PK PEEK anchors are manufactured with polyether ether ketone (PEEK), as per the ASTM F2026 standard.

CordFit #2 (USP standard) sutures are manufactured with ultra-high-molecular-weight polyethylene (UHMWPE). Sutures used may be blue, white and/or white-blue. The Cobalt Chromite Blue-Green Spinel pigment is used for the blue threads. White threads are not pigmented.

The 30" TapeFit Max sutures are manufactured with a mixture of ultra-high-molecular-weight polyethylene (UHMWPE) and polypropylene, for white-blue suture, or a mixture of ultra-high-molecular-weight polyethylene (UHMWPE) and polyamide (nylon), for white-black suture.

The blue pigment used for these sutures is phthalocyanine, while the black pigment used is hematein. White fibers in the suture are not pigmented.

The insertion device is made up of the following materials:

- Handle: it is manufactured with medical grade ABS. There are different colors to identify the anchor model for which the insertion device is used. The colors used to identify the products are mentioned below:
  - Knotix PK is identified with a white and/or cyan blue color,
  - Knotix SPD PK is identified with a white and/or green color,
  - SpinFit PK features a white, orange and/or green handle, and
  - QuickFit PK comes in white, green and/or orange.
- Capillary: manufactured with stainless steel, as per the AISI 304 standard.
- Rod: manufactured with stainless steel, as per the AISI 304 standard.

Check the packaging labels for the anchor materials.

### 13. LEARN MORE

To learn more about this product, visit our website:  
[www.fergusmedical.com](http://www.fergusmedical.com)

### 14. SYMBOLS USED IN LABELS

	<b>CATALOG NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>SERIAL NUMBER</b>
	<b>PRECAUTION</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>DATE OF MANUFACTURE</b>
	<b>STERILE STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>DOUBLE STERILE BARRIER SYSTEM</b>
	<b>MANUFACTURER</b>
	<b>CHECK INSTRUCTIONS FOR USE</b>
	<b>DO NOT USE IF PACKAGING IS DAMAGED</b>
	<b>DO NOT RE-STERILIZE</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>KEEP DRY</b>
	<b>MR-SAFE</b>

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**PROMEDON S.A.**