# Knotix Ti

# Knotix Ti 2.8

---

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ESPAÑOL

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PORTUGUÊS

Instructions for use

ENGLISH

M-351-00204 (Rev. 05) / 10.04.2025

# Fergus

Fergus

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## 1. DESCRIPTION OF THE DEVICE

The Fergus titanium anchors with sutures include the Knotix Ti and Knotix Ti 2.8 models. These are titanium devices that are designed to be threaded and feature eyelets to pass the suture. They are supplied sterile, threaded with non-absorbable sutures and inside a disposable insertion device to aid in the placement of the implant.

The titanium anchors with sutures are designed to attach soft tissue to the bone by tying knots.

The Knotix Ti model features a threaded anchor and an internal eyelet where two CordFit #2 sutures are passed (USP standard). Sutures may come in blue, white or white-blue colors. The four tips of these two sutures are connected to the suture retention area located in the handle of the insertion device. This anchor model is available in 5.5 mm or 6.5 mm diameters with 16.3 mm in length.

Figure 1: Knotix Ti anchor (illustrative picture).  
1: Suture retention area.  
2: Titanium anchor.

The Knotix Ti 2.8 model features a threaded anchor and includes a single CordFit #2 suture (USP standard). The suture may be blue, white or white-blue, and it comes threaded in the anchor eyelet. The ends of these two sutures are connected to the suture retention area located in the handle of the insertion device. This anchor model is available with a diameter of 2.8 mm and a length of 11.6 mm.

Figure 2: Knotix Ti 2.8 anchor (illustrative picture).  
1: Suture retention area.  
2: Titanium anchor.

## 2. HOW SUPPLIED

Contents:

- **Knotix Ti Model:** implantable titanium anchor with two CordFit #2 sutures and an insertion device.
- **Knotix Ti 2.8 Model:** implantable titanium anchor with one CordFit #2 suture and an insertion device.

Find the model, anchor size, materials, and suture and insertion device color on each insertion device label.

All components are supplied sterile and ready-to-use.

The availability of models and indications for use may vary from one market to another, which is subject to regulatory approvals in each country.

Please, check the Indications for use section to find the indications authorized for each country.

If you have any questions, please direct them to your commercial representative.

### 3. INDICATIONS FOR USE

Knotix Ti and Knotix Ti 2.8 titanium suture anchors are intended to provide fixation of soft tissue to bone within the body joints. They are designed to be used in the shoulder, wrist, knee and ankle.

**NOTE:** In Mexico, they are authorized to be used in the shoulder and ankle only.

### 4. CONTRAINDICATIONS

The device shall not be implanted in the following cases:

- Quality or pathologic conditions of the soft tissue to be fixed which could affect the safety of fixation with the use of sutures. The physician should assess bone quality carefully in patients with lack of bone maturity before starting the orthopedic surgical procedure. The growth plate should not be deviated, interrupted or altered.
- Poor quality or pathological conditions of the soft tissue to be repaired, which could impair secure fixation by suture.
- Limited blood irrigation or previous infection which may delay healing.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions of their physical activity during the osseointegration period.
- Surgical interventions not listed in section "Indications for Use."
- Inadequate neuromuscular status.

### 5. WARNINGS

- Before using this product, read the instructions for use fully and thoroughly.
- The device must be used only by a medical professional trained in the surgical techniques indicated for the application of these implants.
- The law restricts the use of this device to physicians or people authorized by a physician.
- DO NOT use the product if the packaging is open or damaged, as product sterilization may have been compromised.
- Products are designed for a SINGLE use. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RSTERILIZED, since this could damage their performance and increase the risk of inadequate sterilization and cross contamination.

- Maximum precautions must be taken to avoid contamination. Operating room conditions must meet hospital, administrative and/or local government requirements.
- Do not use the device after the expiration date stated on the label. The performance, safety and/or sterility of the device cannot be guaranteed after the expiration date.
- Ensure that the product is disposed of following local regulations and be aware of the environmental contamination risks to the staff and the patients.
- The product must be stored in its original packaging.
- The patient's sensitivity to device materials must be considered before implantation.
- All surgical instruments used in the intervention meant for the application of this implant must be sterile and fit for use in aseptic environments.
- Product presentation should not be altered without following the surgical technique, as this could impair the implant's performance.
- If the implant disengages from the insertion device, **DO NOT REASSEMBLE**, since this could damage its performance and increase the risk of inadequate sterilization.
- The fixation provided by the implant must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing the product and the bone.
- Patients will receive detailed instructions about the use and limitations of the product.
- If the decision is made to explant the product, the risks of undergoing a second surgical intervention should be considered. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.
- Open devices and implants should be discarded, even if they have not been used or handled, since this could increase the risk of inadequate sterilization and cross contamination.
- Biohazardous waste such as explanted devices, sutures, needles and contaminated surgical material must be disposed of safely and following the policy of the medical facility.
- Do not use the product if any faults are detected.
- Do not use products or sutures that are not recommended by the manufacturer, since implant performance for the repair cannot be guaranteed in this case.
- The manufacturer is not responsible for failures or damage caused when these devices are combined with products not recommended by the manufacturer.
- The manufacturer shall not be liable for any effect or consequence that may result from not following specific surgical techniques or instructions, such as a mistaken indication, an inadequate use of a surgical technique, asepsis issues, etc., which could affect patient safety and device performance.

## 6. PRECAUTIONS

- The adequate implant model, quantity and size should be determined at the surgeon's professional discretion based on the specific indication and the patient's medical history.
- Surgeons are recommended to check the specific surgical technique for the product before proceeding with the operation.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper device selection and installation, are essential for satisfactory use of the implant.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.
- The surgeon must consider the manufacturer-recommended installation systems in order to implant the product properly.
- Do not remove the suture until the anchor has been completely inserted into the desired site, since this could impair implantation.
- Using excessive force during insertion may cause faults in the suture anchor or the insertion device.
- Sutures can break if handled or managed with sharp or cutting instruments.
- The decision of removing the implant shall be made considering the potential risk of a second surgical intervention for the patient. After the explantation, appropriate post-operative controls must be carried out.

## 7. SUPPLY AND STORAGE

- Titanium anchors with sutures are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- In order to maintain product quality and effectiveness, it is marketed in a double Tyvek pouch (primary packaging). After sterilization, the product is protected during transport and storage by placing it inside a cardboard box (external packaging). Additionally, the product is supplied in a plastic protective wrap covering the implant, which is used to prevent Tyvek pouch damage. The product is identified with labels.
- Keep the product away from any moisture or heat sources.
- The product is recommended to be kept at room temperature.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- Product shelf life is 3 years.

## **8. MAGNETIC RESONANCE (MR) ENVIRONMENT**

MR conditional.

Knotix Ti and Knotix Ti 2.8 titanium sutures anchors by Fergus, made of Ti-6Al-4V ELI titanium alloy as per ASTM F136, feature a low magnetic susceptibility and electric conductivity. Therefore, there is a low risk of displacements due to induced magnetic forces or torsion, and of heating induced by radiofrequency or imaging devices are low.

Patients who have been implanted with these products may be subjected to magnetic resonance screenings safely under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, only.
- Magnetic field with a maximum spatial gradient of 3000 Gauss/cm.
- Maximum reported by the Magnetic Resonance System, specific absorption rate (SAR) averaged all over the body at 1 W/kg for 15 minutes of exploration in the normal mode of the MR system.

Image quality could be compromised in the function area or its surroundings. Magnetic resonance parameters need to be optimized to compensate for the presence of the implant.

## **9. POTENTIAL COMPLICATIONS**

Any surgical intervention involves a certain level of risk. Some complications that may occur are:

- Deep or superficial infections
- Atypical reactions to foreign bodies
- Pain at the incision or intervention site
- Implant loosening or detachment
- Bone injuries
- Incomplete or inadequate lesion repair or lesion recurrence

## **10. INFORMATION FOR PATIENTS**

Surgeons are responsible for informing patients or their representatives, prior to the surgery, of the potential complications associated to the device implantation (see section "Potential complications").

The company and their distributors delegate to the surgeon the responsibility to inform the patient of the advantages and possible risks related to the implantation and use of the product.

Postoperative care is important. The patient will be instructed by the physician on the limitations of the implant and will be warned about the care needed for healing.

Patients should be warned that if they decide to explant the product, they should consider the risks of undergoing a second surgical procedure. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.

The performance of additional surgery does not guarantee the total repair of complications.

11. SURGICAL TECHNIQUE

The surgeon must select the appropriate approach based on the repair to be performed. The surgeon may opt for an open, minimally-invasive, or arthroscopic approach, as preferred by the surgeon.

Additionally, the surgeon must choose the number of implants and the appropriate fixation layout, depending on the lesion type, the anatomical region, the patient’s medical history and other criteria deemed to be important for the successful use of the product in the surgical intervention.

The following instructions for surgical techniques only refer to the device implantation; they do not include other procedures that are essential to perform a surgical intervention, such as patient preparation, access portals, surgery completion, among others.

Some interventions associated to each titanium anchor model are suggested below. However, other procedures may be performed if the physician deems it appropriate as long as the indications for use are followed:

***Knotix Ti***

Knotix Ti is recommended for:

|          |                                                                                                                                                                                                                                                         |
|----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Shoulder | <ul style="list-style-type: none"><li>• Rotator cuff repair</li><li>• Proximal biceps tenodesis</li></ul>                                                                                                                                               |
| Knee     | <ul style="list-style-type: none"><li>• Lateral extra-articular tenodesis (LET)</li><li>• Posterolateral complex reconstruction</li><li>• Posteromedial complex reconstruction</li><li>• Medial patellofemoral ligament (MPFL) reconstruction</li></ul> |
| Ankle    | <ul style="list-style-type: none"><li>• Lateral instability repair</li><li>• Medial instability repair</li></ul>                                                                                                                                        |

Steps for the right implantation of the Knotix Ti anchor:

- 1)** Position the anchor in the desired implantation site.
- 2)** Rotate the (insertion device) handle clockwise to enable anchor implantation in bone tissue by threading. While turning the handle, a slight forward pressure (in the proximal and distal direction) is recommended to make implantation easier.

ATTENTION: It is important to ensure that the anchor has been fully inserted to avoid possible damages to surrounding tissues and/or implant detachments. To this end, threading should continue until the laser line located in the distal end of the insertion device shaft has fully entered the hole where the anchor is placed.

- 3)** Release sutures from the handle retention area and remove the insertion device.

**4)** Use a suture passer to pass sutures through the tissue, tie them and cut off the excess suture.

**ATTENTION:** If an open approach is used, the suture must be threaded through a needle to pass through the tissue, which allows for tissue fixation to the anchoring point. Tie surgical knots to end the fixation.

**NOTE:** The needle is not included in the kit.

After use, dispose of the product, sutures, needles, contaminated disposable surgical material and packaging, in accordance with hospital, administrative, or local government procedures.

## ***Knotix Ti 2.8***

Knotix Ti 2.8 is recommended for:

|          |                                                                                                                  |
|----------|------------------------------------------------------------------------------------------------------------------|
| Shoulder | <ul style="list-style-type: none"><li>• Glenoid labral repair</li></ul>                                          |
| Wrist    | <ul style="list-style-type: none"><li>• Triangular ligament reconstruction or repair</li></ul>                   |
| Ankle    | <ul style="list-style-type: none"><li>• Lateral instability repair</li><li>• Medial instability repair</li></ul> |

Steps for the right implantation of the Knotix Ti 2.8 anchor:

**1)** Identify the anchor insertion site. For shoulder instability repair only, place the 3.2-mm implantation guide (FG-IN013000) in the insertion site. Then, insert the anchor through the guide.

**NOTE:** If using an open approach, the use of the guide is not required

**2)** Rotate the (insertion device) handle clockwise to enable anchor implantation in bone tissue by threading. While turning the handle, a slight forward pressure (in the proximal and distal direction) is recommended to make implantation easier.

**ATTENTION:** It is important to ensure that the anchor has been fully inserted to avoid possible damages to surrounding tissues and/or implant detachments. To this end, threading should continue until the laser line located in the distal end of the insertion device shaft has fully entered the hole where the anchor is placed.

**3)** Release sutures from the handle retention area and remove the insertion device and then the implantation guide.

**4)** Use a suture passer to pass sutures through the tissue, tie them and cut off the excess suture.

**ATTENTION:** If an open approach is used, the suture must be threaded through a needle to pass through the tissue, which allows for tissue fixation to the anchoring point.

**5)** Tie surgical knots to end the fixation.

**NOTE:** The needle is not included in the kit.

After use, dispose of the product, sutures, needles, contaminated disposable surgical material and packaging, in accordance with hospital, administrative, or local government procedures.

## 12. MATERIAL SPECIFICATIONS

Knotix Ti and Knotix Ti 2.8 titanium anchors with sutures are manufactured using a Ti-6Al-4V ELI titanium alloy as per the ASTM F136 standard.

CordFit #2 sutures (USP standard) are made of ultra-high-molecular-weight polyethylene (UHMWPE). Sutures can be blue, white or white/blue in color. The *Cobalt Chromite Blue-Green Spinel* pigment is used for the blue threads. White threads are not pigmented.

The insertion device consists of two parts:

- **Handle:** It is made of medical grade ABS and it has different colors to identify the anchor model for which this insertion device is used. A white and/or cyan blue device is used for a Knotix Ti suture anchor, whereas a white and/or green insertion device is used for a Knotix Ti 2.8 suture anchor.
- **Shaft:** It is made of stainless steel as per AISI 304.

See the packaging label to learn more about the anchor materials.

## 13. MORE INFORMATION

If you need more information about this product, visit our web page, [www.fergusmedical.com](http://www.fergusmedical.com)

## 14. SYMBOLS USED IN THE LABELS

**CATALOGUE NUMBER**

**BATCH NUMBER**

**SERIAL NUMBER**

**PRECAUTION**

**EXPIRATION DATE**

**DO NOT REUSE**

**DATE OF MANUFACTURE**

**STERILE. STERILIZATION METHOD:  
ETHYLENE OXIDE**

**DOUBLE STERILE BARRIER SYSTEM**

**MANUFACTURER**

**SEE INSTRUCTIONS FOR USE**

**DO NOT USE IF THE PACKAGING IS  
DAMAGED**

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**PROMEDON S.A.**