# Suturas especiales

---

**CordFit** Infinity

**CordFit** Stiff

**CordFit** with Tapered Needle

**CordFit** Meniscal Repair

**CordFit** A Infinity

**CordFit** A Stiff

**CordFit** A with Tapered Needle

**CordFit** A Meniscal Repair

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ESPAÑOL

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PORTUGUÊS

Instructions for use

ENGLISH

M-350-00735 (Rev. 04) / 26.11.2024

Fergus

Fergus

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**RESSONANCIA MAGNÉTICA SEGURA**

## 1. PRODUCT DESCRIPTION

The Fergus-brand special suture line is made up of the CordFit Infinity, CordFit Stiff, CordFit with Tapered Needle, CordFit Meniscal Repair, CordFit A Infinity, CordFit A Stiff, CordFit A with Tapered Needle and CordFit A Meniscal Repair. These are surgical, sterile, non-absorbable sutures manufactured with braided multifilaments made of ultra-high-molecular-weight polyethylene (UHMWPE) or a mixture of UHMWPE and polypropylene or polyamide which may be attached posteriorly to the needle, depending on the model. The length of the suture thread is 38" (965 mm) for all special suture models.

The CordFit Infinity, CordFit A Infinity, CordFit with Tapered Needle, CordFit A with Tapered Needle, CordFit Meniscal Repair and CordFit A Meniscal Repair models are designed for applications involving soft tissue attachment and/or approximation.

The CordFit Meniscal Repair and CordFit A Meniscal Repair models are specifically designed for repairing the meniscus.

The CordFit Stiff and CordFit A Stiff models are designed to assist in procedures, mainly for the passage of suture through bone tunnels and suture recovery out of the joint area.

The CordFit Stiff and CordFit A Stiff special sutures are the only models that have not been designed for implantation and are for temporary use exclusively.

The CordFit Infinity, CordFit Stiff, CordFit with Tapered Needle and CordFit Meniscal Repair models are thread-type UHMWPE sutures that come in white, white-blue and blue colors.

- #2 UHMWPE CordFit Infinity suture, white, looped, with straight needle and nitinol loop.
- #2 UHMWPE CordFit Infinity suture, white-blue, looped, with straight needle and Nitinol loop.
- #2 UHMWPE CordFit Infinity suture, blue, looped, with straight needle and Nitinol loop.
- #2 UHMWPE CordFit Infinity suture, white, looped, with straight needle and stainless steel loop.
- #2 UHMWPE CordFit Infinity suture, white-blue, looped, with straight needle and stainless steel loop.
- #2 UHMWPE CordFit Infinity suture, blue, looped, with straight needle and stainless steel loop.

- #2 UHMWPE CordFit with Tapered Needle suture, white, with 1/2 circle needle.
- #2 UHMWPE CordFit with Tapered Needle suture, white-blue, with 1/2 circle needle.
- #2 UHMWPE CordFit with Tapered Needle suture, blue, with 1/2 circle needle.
- #2 UHMWPE CordFit Stiff suture, white, with stiffened end.
- #2 UHMWPE CordFit Stiff suture, white-blue, with stiffened end.
- #2 UHMWPE CordFit Stiff suture, blue, with stiffened end.
- #0 UHMWPE CordFit Meniscal Repair suture, white-blue, with two needles.
- #0 UHMWPE CordFit Meniscal Repair suture, blue, with two needles.
- #2-0 UHMWPE CordFit Meniscal Repair suture, white-blue, with two needles.
- #2-0 UHMWPE CordFit Meniscal Repair suture, blue, with two needles.

This suture line includes totally- or partially-pigmented threads with a specific pattern. The blue color of the suture threads is achieved by using the Cobalt Chromite Blue-Green Spinel (C.I. Pigment Blue 36) pigment. White threads are not pigmented.

The CordFit A Infinity, CordFit A Stiff, CordFit A with Tapered Needle and CordFit A Meniscal Repair models are thread-type sutures manufactured using UHMWPE and polypropylene (white-blue) or UHMWPE and polyamide (white-black).

- #2 UHMWPE and polypropylene CordFit A Infinity suture, white-blue, looped, with straight needle and Nitinol loop.
- #2 UHMWPE and polyamide CordFit A Infinity suture, white-black, looped, with straight needle and Nitinol loop.
- #2 UHMWPE and polypropylene CordFit A Infinity suture, white-blue, looped, with straight needle and stainless steel loop.
- #2 UHMWPE and polyamide CordFit A Infinity suture, white-black, looped, with straight needle and stainless steel loop.
- #2 UHMWPE and polypropylene CordFit A with Tapered Needle suture, white-blue, with 1/2 circle needle.
- #2 UHMWPE and polyamide CordFit A with Tapered Needle suture, white-black, with 1/2 circle needle.
- #2 UHMWPE and polypropylene CordFit A Stiff suture, white-blue, with stiffened end.

- #2 UHMWPE and polyamide CordFit A Stiff suture, white-black, with stiffened end.
- #2-0 UHMWPE and polypropylene CordFit A Meniscal Repair suture, white-blue, with two needles.
- #2-0 UHMWPE and polyamide CordFit A Meniscal Repair suture, white-black, with two needles.

This suture line includes partially-pigmented threads with a specific pattern. In white-blue sutures, the blue color is achieved with the phthalocyanine pigment, while in the case of white-black sutures the black color is achieved with the hematein pigment. White threads are not pigmented.

## 2. PRODUCT PRESENTATION

Special sutures from the Fergus line are available in the following presentations:

- CordFit with Tapered Needle y CordFit A with Tapered Needle sutures: supplied in 6- or 12-unit boxes. Each product unit is packaged individually.
- CordFit Infinity y CordFit A Infinity sutures: supplied in 6- or 12-unit boxes. Each product unit is packaged individually.
- CordFit Stiff y CordFit A Stiff sutures: supplied in 3- or 6-unit boxes. Each product unit is packaged individually.
- CordFit Meniscal Repair y CordFit A Meniscal Repair sutures: supplied in 3- or 6-unit boxes. Each product unit is packaged individually.

Check the individual product labels for the model, material, color and dimensions. All components are provided sterile and ready for use.

## 3. INDICATIONS FOR USE

Fergus-brand sutures are intended for applications that involve approximating or joining soft tissue, and can be used as components in procedures that include repairs using allograft or autograft tissue in arthroscopic surgery, minimally-invasive orthopedic surgery and/or open orthopedic surgery.

CordFit Stiff and CordFit A Stiff are specifically indicated as temporary components to assist in procedures involving repairs using allografts or autografts, arthroscopic surgery, minimally-invasive orthopedic surgery, and/or open orthopedic surgery.

#### **4. CONTRAINDICATIONS**

- The products must not be used for heart surgery, wound closing or surgical procedures not mentioned in the “Indications for Use” section.
- The sutures are not to be implanted in the following cases (not applicable for CordFit Stiff and CordFit A Stiff models, which are not indicated to be implanted in any case):
  - Tissue quality or pathologies of the tissue to be reconstructed which may affect the safety of the suture attachment.
  - Limited blood irrigation and/or previous infection which may delay healing and/or consolidation.

#### **5. WARNINGS**

- Read the instructions for use carefully and completely before use.
- The product must be used only by a medical professional trained in the surgical techniques indicated for the use of these products.
- The law allows this product to be used only by physicians or people authorized by a physician.
- **DO NOT** use the product if the packaging is open or damaged, since it may not be possible to guarantee product safety and effectiveness.
- Do not use the product if any faults are detected.
- The product has been designed for **SINGLE** use only. Therefore, **IT MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED**, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative or governmental requirements.
- Do not use the product after the expiration date that appears on the label. Product performance, safety and/or sterility cannot be guaranteed after this date has passed.
- Do not expose the product to a direct heat source.
- Make sure that the product is disposed of as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- The product must be stored in the original packaging.
- All surgical instruments used in the surgical procedure to install this product must be sterile and suitable for use in aseptic environments.
- Once the product packaging is open, all sutures that do not remain implanted in the patient must be discarded, regardless of them not being used or manipulated.

Additionally, the protective tube and the suture of the CordFit Stiff and CordFit A Stiff models, as well as the needles that have not been used, must be properly disposed of to reduce the risk of improper sterilization and cross-contamination.

- Bio-hazardous waste, such as explanted sutures, needles, and contaminated surgical materials, must be disposed of safely and as per site policy.
- The suture has not been found or proven to be safe and effective for use as an artificial tendon or ligament.
- The manufacturer shall not be held liable for faults or damages caused when these products are combined with other products not recommended by the manufacturer.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, etc., which may affect patient safety and product performance.

The following warnings are not applicable for CordFit Stiff and CordFit A Stiff models, which are not indicated to be implanted in any case:

- Patient sensitivity to product materials must be considered before implantation.
- The fixation provided by the suture must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing of the suture.
- If the suture is to be explanted, the risks of the possibility of a second procedure on the patient must be considered. Suture extraction must be accompanied by an appropriate post-operative approach.
- Needles should be handled carefully to avoid accidental punctures.

## **6. PRECAUTIONS**

- Surgeons must follow their professional judgment to determine the appropriate product model based on the specific indication and the history of the patient.
- Surgeons are advised to review the product-specific surgical technique before the procedure.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper product selection, are essential for satisfactory use of the product.

- Sutures may break if handled or controlled using sharp, pointed or cutting instruments.

The following precautions are not applicable for CordFit Stiff and CordFit A Stiff models, which are not indicated to be implanted in any case:

- Do not hold the needles by their ends with forceps or pliers to avoid damaging them.
- Do not warp needles or use needles that are warped, as they may lose their strength or break more easily.
- Avoid damage to adjacent tissues or accidental punctures due to improper needle handling.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.

## **7. SUPPLY AND STORAGE**

- The product is supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- The suture can be supplied with a needle, two needles, or without a needle, depending on the model.
- In order to maintain product quality and effectiveness, the CordFit Infinity, CordFit A Infinity, CordFit with Tapered Needle, CordFit A with Tapered Needle, CordFit Meniscal Repair and CordFit A Meniscal Repair models are supplied inside a protection card and in a double Tyvek pouch. CordFit Stiff and CordFit A Stiff models are supplied in a double Tyvek pouch. CordFit Stiff and CordFit A Stiff models are supplied in a double Tyvek pouch. After the sterilization, the product is protected during transportation and storage by placing it inside a car dboard box. Labels may be added for identification.
- Keep the product away from any sources of heat or humidity.
- Room temperature storage is advised.
- **IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.**
- **DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.**
- Product shelf life is 3 years.

## **8. MAGNETIC RESONANCE (MR) ENVIRONMENTS**

- **Conditional MR Compatibility and MR Safety**  
Sutures are made of ultra-high-molecular-weight polyethylene (UHMWPE) or a mixture of UHMWPE with polypropylene or polyamide, which makes them MR safe.

## 9. POSSIBLE COMPLICATIONS

All surgical procedures involve a certain risk. Some of the complications that could arise are:

- Deep and superficial infections.
- Possible allergies and other reactions to device materials.
- Pain on the site of the incision or the procedure.
- Edema and erythema on the site of the wound.
- Tissue lesions.
- Incomplete or inappropriate lesion correction, or lesion recurrence

The following complications are not applicable for CordFit Stiff and CordFit A Stiff models, which are not indicated to be implanted in any case:

- Reactions to foreign bodies
- Wound dehiscence
- Gradual suture encapsulation in fibrous connective tissue

## 10. PATIENT INFORMATION

Surgeons are responsible for informing patients and their representatives before the surgery about the possible complications related to the use or implantation of the product (see the "Possible Complications" section).

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated with the implantation and use of the product.

Post-operative care is important after suture implantation (this is not applicable for CordFit Stiff and CordFit A Stiff models, since they are not implantable). Patients will receive instructions from the physician about the limitations required by the implant and will be warned about the necessary precautions they must apply for healing.

The patient must be warned that if the decision is made to explant the product (except for the CordFit Stiff and CordFit A Stiff, which are not implantable), the risks of undergoing a second surgical procedure should be considered. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach. Performing additional surgeries does not guarantee a complete repair of the complications.

## 11. SURGICAL TECHNIQUE

The surgeon must be trained on the use of accepted techniques for surgical knots, number of sutures and

configurations, on the surgical interventions mentioned in the section “Indications for use”, the type of injury to repair, the clinical history of the patient and other criteria that the professional considers important to achieve a satisfactory use of the product in the surgical intervention.

Some interventions associated to each special sutures model are suggested below. However, other procedures may be performed if the physician deems it appropriate as long as the indications for use are followed:

- CordFit with Tapered Needle and CordFit A with Tapered Needle.

They are recommended for the suture of ligaments, tendons, or grafts in orthopedic surgery.

- CordFit Infinity and CordFit A Infinity.

They are recommended for the preparation of grafts in orthopedic surgery.

- CordFit Stiff and CordFit A Stiff.

They are recommended to be used as accessories for supporting procedures such as the reconstruction of the anterior and/or posterior cruciate ligament, either by means of an all-inside or retroconstruction technique, through a bone tunnel. They may also be used to repair the meniscal root. None of these indications involve implanting the product.

- CordFit Meniscal Repair and CordFit A Meniscal Repair.

They are recommended to repair the meniscus by means of the inside-out technique.

This product could be used in combination with other medical products, but the surgical technique associated to those products must be followed.

After use, discard the product and its packaging per local hospital, administrative or governmental procedures.

## 12. MATERIAL SPECIFICATIONS

The special sutures may be manufactured with:

- Ultra-high-molecular-weight polyethylene (UHMWPE): sutures may be white, white-blue, or blue. The blue color of the suture threads is achieved by using the Cobalt Chromite Blue-Green Spinel pigment. White threads are not pigmented.
- UHMWPE and polypropylene: sutures may be white-blue. The pigment used for the blue threads of the suture is phthalocyanine. White fibers in the suture are not pigmented.
- UHMWPE and polyamide: sutures may be white-black. The pigment used for the black threads of the suture is hematein. White fibers in the suture are not pigmented.

The needles are manufactured with AISI 300 series stainless steel.

The protective tube for the CordFit Stiff and CordFit A Stiff sutures is made of polypropylene.

CordFit Infinity and CordFit A Infinity sutures feature a Nitinol or stainless steel loop crimped to the needle.

Check the packaging labels for more information.

### 13. LEARN MORE

To learn more about this product, visit our website: [www.fergusmedical.com](http://www.fergusmedical.com)

### 14. SYMBOLS USED IN LABELS

	<b>CATALOG NUMBER</b>
	<b>BATCH NUMBER</b>
	<b>SERIAL NUMBER</b>
	<b>PRECAUTION</b>
	<b>EXPIRATION DATE</b>
	<b>DO NOT REUSE</b>
	<b>MANUFACTURING DATE</b>
	<b>STERILE. STERILIZATION METHOD: ETHYLENE OXIDE</b>
	<b>DOUBLE STERILE BARRIER SYSTEM</b>
	<b>MANUFACTURER</b>
	<b>CHECK INSTRUCTIONS FOR USE</b>
	<b>DO NOT USE IF PACKAGING IS DAMAGED</b>
	<b>DO NOT RE-STERILIZE</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>KEEP DRY</b>
	<b>MR-SAFE</b>

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**PROMEDON S.A.**