Instrumentales estériles

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Instructions for use ENGLISH

M-350-00622 (Rev. 03) / 28.06.2024

Fergus

ESPAÑOL

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ENGLISH

1. DESCRIPTION OF THE DEVICE

The devices that conform the Fergus sterile instrument system are single-use sterile manual orthopedic surgical instruments.

These devices may include instruments used for bone perforation, soft tissue passage and suture retrieval. The instruments have an ergonomic, low-profile design intended to gain manual control from different positions and for different functions.

This system includes instruments to perform bone perforations, such as the 4 mm spade-tip pin with open eyelet and the 4 mm spade-tip pin with closed eyelet, to pass through soft tissue, such as the needle for SelfPass forceps, and to retrieve suture, such as the Nitinol loop.

The instruments are designed to simplify surgical procedures for shoulder, knee, foot and ankle reconstruction.

2. INDICATIONS FOR USE

The Fergus sterile instruments are indicated for shoulder, knee, foot and ankle arthroscopic surgical procedures, minimally invasive orthopedic surgery and/or open orthopedic surgery.

3. CONTRAINDICATIONS

No contraindications were found for these devices.

4. WARNINGS

- Read the instructions for use carefully and completely before use.
- The device must be used only by a medical professional trained in the surgical techniques indicated.
- The law restricts the use of this device to physicians or people authorized by a physician exclusively.
- These sterile instruments have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- DO NOT use the device if the packaging is open or damaged, since product sterility will be lost.

- The device must not be used for surgical procedures that are not included in the "Indications for Use" section.

- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date.

- The product must be stored in the original packaging.

- To prevent damaging the Needle for SelfPass forceps, do not use it outside the internal joint area.

- Avoid damaging neurovascular structures when handling these devices.

- None of these devices should be used by means of impact.

- To prevent damage or breakage during use, do not apply excessive force on the spade-tip pin. Do not use the spade-tip pin as a lever.

- Do not implant. These devices are not designed for implantation.

- The Spade-tip pin and the Needle for SelfPass forceps are extremely sharp. Handle these items with care to prevent possible injury.

- The spade-tip pin must be removed before flexing the joint, since the device could bend or break.

- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be according to local hospital, administrative and/or governmental requirements.

- Make sure that the product is disposed as per local regulations, and be wary of the risks of environmental, staff and patient contamination.

- Bio-hazardous waste, such as explanted devices, sutures, instruments and contaminated surgical materials, must be disposed of safely and as per site policy.

- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the manufacturer.

5. PRECAUTIONS

- Before using the product, remove the packaging and the tip protector, if used.

- After the instrument has been inserted in the joint, do not apply excessive force or perform any

inappropriate movements, since it could break and become lodged in the soft tissue or disappear from the arthroscopic view of the surgical field.

- Before using the device, examine it for damage. Do not use if damaged.

6. SUPPLY AND STORAGE

- The instruments are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.

- In order to maintain product quality and effectiveness, it is marketed in a double Tyvek pouch (primary packaging). After sterilization, the product is protected during transport and storage by placing it inside a cardboard box (external packaging). Additionally, the product may be supplied in a protective wrap covering the implant, which is used to prevent Tyvek pouch damage. The product is identified with labels.

- These devices must be stored in their original unopened packaging, and they must not be used after their expiration date.

- Keep the product away from any sources of heat or humidity.

- Room temperature storage is advised for this product.






- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.

- Product shelf life is 3 years.

7. POSSIBLE COMPLICATIONS

No possible complications were found for these devices.

8. SYMBOLS USED IN LABELS

	CATALOG NUMBER
	BATCH NUMBER
	SERIAL NUMBER
	PRECAUTIONS
	EXPIRATION DATE
	DO NOT REUSE
	MANUFACTURING DATE
	STERILE. STERILIZATION METHOD: ETHYLENE OXIDE
	DOUBLE STERILE BARRIER SYSTEM
	MANUFACTURER
	DO NOT USE IF PACKAGING IS DAMAGED
	CHECK INSTRUCTIONS FOR USE
	DO NOT RE-STERILIZE
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY

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PROMEDON S.A.