Spirom PK
Tenodesis

SpinFit PK
Tenodesis

Anclajes de PEEK para tenodesis

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Instructions for use
ENGLISH

M-350-00192 (Rev. 02) / 18.12.2024

Fergus

Fergus

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MANTER SECO

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1. DEVICE DESCRIPTION

The Fergus PEEK tenodesis anchors, Spirom PK Tenodesis and SpinFit PK Tenodesis models, are polyetheretherketone (PEEK) implants supplied sterile and they come with or without a disposable insertion device.

The Fergus PEEK tenodesis anchors are designed for soft tissue fixation to bone.

The Spirom PK Tenodesis is a hexagonal-cored screw manufactured in PEEK and it is available in the following dimensions:

- 6 mm diameter and 19 mm length.
- 7 mm diameter and 19 mm length.
- 8 mm diameter and 19 mm length.

The dimensions are shown in the product labels.

The implantation must be performed with a reusable hex screwdriver.


Figura 1: Spirom PK Tenodesis.

SpinFit PK Tenodesis consists of a screw and an eyelet used to fix soft tissue inside the hole in the bone. These are provided already coupled with a disposable insertion device for their placement.

On its distal end, SpinFit PK Tenodesis presents two types of eyelets:

- I. Suture eyelet: white-blue, cord-type #2 (USP standard) ultra-high-molecular-weight polyethylene (UHMWPE) suture eyelet, which is inserted inside the capillary.
- II. Fork eyelet: bifurcated PEEK eyelet which is detached from the threaded implant and attached to the insertion device through a white-blue, cord-type #2 (USP standard) UHMWPE retention suture.

This product is available in the following dimensions:

- 6 mm diameter and 19 mm length (suture eyelet) or 23 mm length (fork eyelet).
- 7 mm diameter and 19 mm length (suture eyelet) or 23 mm length (fork eyelet).
- 8 mm diameter and 19 mm length (suture eyelet) or 23 mm length (fork eyelet).

The dimensions are shown in the product labels.

NOTE. The length stated for the SpinFit PK Tenodesis with fork eyelet equals the addition of the lengths of the screw and the PEEK eyelet.

The soft tissue (graft) is positioned in the hole in the bone using the eyelet. Threading the insertion device makes the implant advance and fixate the soft tissue by means of interference between the bone and the implant.

Figure 2: SpinFit PK Tenodesis (suture eyelet).

Figure 3: SpinFit PK Tenodesis (fork eyelet).

2. PRODUCT PRESENTATION

Contents:

Spirom PK Tenodesis Model:

- One PEEK screw

SpinFit PK Tenodesis Model:

- One PEEK screw with suture eyelet and a disposable insertion device; or
- One PEEK screw with fork eyelet and a disposable insertion device.

Check the individual product labels to find the model, implant size, materials and the eyelet type.

All components are provided sterile and ready for use.

3. INDICATIONS FOR USE

The Fergus PEEK tenodesis anchors are indicated for the fixation of soft tissue, including ligaments and tendons, to bone. They are designed to be used in the shoulder, knee and ankle.

4. CONTRAINDICATIONS

The device should not be implanted in the following cases:

- Insufficient bone quality or quantity for attachment. The physician must carefully evaluate the bone quality of patients with bone immaturity before starting with the orthopedic surgical procedure, and must avoid deviating, interrupting or altering the growth cartilage.
- Quality or pathologic conditions of the soft tissue to be fixed which could affect the safety of the fixation with the implant.
- Limited blood irrigation or previous infection which may delay healing.
- Conditions that may limit the ability and/or the will of the patient to follow the physicians' indications regarding post-operative care, which may include significant restrictions of their physical activity.
- Bone-to-bone fixation and other surgical procedures not mentioned in the "Indications for Use" section.
- Inappropriate neuromuscular conditions.

5. WARNINGS

- Read the instructions for use carefully and completely before use.
- The device must be used only by a medical professional trained in the surgical techniques indicated for placing these implants.
- The law restricts the use of this device to physicians or people authorized by a physician.
- DO NOT use the product if the packaging is open or damaged since product sterility may be lost.
- The products have been designed for SINGLE use only. Therefore, THEY MUST NOT BE REUSED, REPROCESSED OR RESTERILIZED, since this could impair their effectiveness and increase the risk of inappropriate sterilization and cross-contamination.
- The strictest precautions must be implemented to prevent contamination. Operating room conditions must be compliant with local hospital, administrative and/or governmental requirements.

- Do not use the device after the expiration date that appears in the label. Device performance, safety and/or sterility cannot be guaranteed after the expiration date has passed.
- Make sure that the product is disposed of as per local regulations, and be wary of the risks of environmental, staff and patient contamination.
- The product must be stored in the original packaging.
- Patient sensitivity to device materials must be considered before implantation.
- All surgical instruments used in the procedure to install this implant must be sterile and suitable for use in aseptic environments.
- Product presentation must not be altered in a way different from what is stated in the surgical technique, since this could affect implant performance.
- For devices with a disposable insertion device, if the implant detaches from the insertion device, DO NOT REATTACH, since this could impair the effectiveness and increase the risk of sterility loss.
- The fixation provided by the implant must be protected after the procedure and until healing is complete. Post-operative guidelines prescribed by the physician must be followed carefully in order to prevent overstressing the product and the bone.
- Patients will receive detailed instructions about the use and limitations of the product.
- If the decision is made to explant the product, the risks of undergoing a second surgical intervention should be considered. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.
- Opened devices and implants must be disposed of, even if they have not been used or handled, since this could increase the risk of inappropriate sterilization and cross-contamination.
- Bio-hazardous waste, such as explants, sutures, instruments and contaminated surgical materials, must be disposed of safely and as per site policy.
- The implant and/or the insertion device may break when the insertion site is not properly prepared with the appropriate instrument before using the implant.
- Do not use the product if any faults are detected.
- Do not use products or sutures that are not recommended by the manufacturer, since implant performance for the repair cannot be guaranteed in this case.

- The manufacturer shall not be held responsible for faults or damages caused when these devices are combined with products not recommended by the manufacturer.
- The manufacturer shall not be held responsible for any effects or consequences that may arise from deviating from these specific techniques or instructions, which may include, but not limited to, an incorrect indication or a misuse of the surgical technique, problems with asepsis, which may affect patient safety and implant performance.

6. PRECAUTIONS

- Surgeons must follow their professional judgment to determine the appropriate implant model, quantity and size based on the specific indication and the history of the patient.
- Surgeons are advised to review the product-specific surgical technique before the procedure.
- Pre-operative and post-operative procedures, including the knowledge about surgical techniques and proper device selection and installation, are essential for satisfactory use of the implant.
- When material sensitivity is suspected, the corresponding tests must be carried out and sensitivity must be ruled out before implantation.
- Product implantation requires preparation in the insertion site.
- The surgeon must consider the manufacturer-recommended installation systems in order to implant the product properly.
- The surgeon must consider the implant size to select the drill bit or bone drill, since an incomplete insertion of the implant may cause a deficient performance of the anchor or an inappropriate implantation.
- If the bone drill or drill bit is not inserted to the appropriate depth, the product may become damaged or broken.
- In order to prevent damage on the implant and/or insertion device, the surgeon must insert it in the same orientation as the hole in the bone.
- Using excessive force during insertion may cause faults in the implant and/or device.
- Insert the capillary of the insertion device into the hole in the bone up to the laser marking of the capillary. Verify and adjust graft tension at this point, since once the implant is placed, the tension cannot be modified.

- Sutures may break if handled or controlled using sharp or cutting instruments.
 - The decision of removing the implant must be made considering the potential risk involved in a second surgical procedure for the patient.
- After the explantation, appropriate post-operative controls must be carried out.

7. SUPPLY AND STORAGE

- PEEK tenodesis anchors are supplied sterile and pyrogen-free. The sterilization method used is ethylene oxide.
- In order to maintain product quality and effectiveness, the product is marketed in a double Tyvek pouch (primary packaging). SpinFit PK Tenodesis features an additional plastic protection covering the tip to prevent damaging the Tyvek pouches. After sterilization, the packed product is protected from light and heat during transport and storage by placing it inside a cardboard box (external packaging). The product is identified with labels.
- Keep the product away from any sources of heat or humidity.
- Room temperature storage is advised.
- IF THE LABEL AND/OR THE PACKAGING ARE OPEN OR ALTERED, DO NOT USE THE PRODUCT.
- DO NOT USE THE PRODUCT AFTER THE EXPIRATION DATE PRINTED IN THE PACKAGING.
- Product shelf life is 3 years.

8. MAGNETIC RESONANCE (MR) ENVIRONMENTS

MR Safety and Compatibility

Tenodesis anchors are manufactured with polyetheretherketone (PEEK); therefore, they are safe for MR.

Sutures are made of ultra-high-molecular-weight polyethylene (UHMWPE) which makes them MR safe.

9. POSSIBLE COMPLICATIONS

All surgical procedures involve a certain risk. Some of the complications that could arise are:

- Deep and superficial infections
- Reactions to foreign bodies
- Pain on the site of the incision or the procedure
- Implant loosening and/or detachment
- Bone lesions
- Incomplete or inappropriate lesion correction, or lesion recurrence.

10. PATIENT INFORMATION

Surgeons are responsible for informing patients and their representatives before the surgery about the possible complications related to the implantation of the product (see the “Possible Complications” section).

The manufacturer and its distributors delegate to the surgeon the responsibility of informing the patient about the advantages and possible risks associated with the implantation and use of the product.

Post-operative care is important. Patients will receive instructions from the physician about the limitations required by the implant and will be warned about the necessary precautions they must apply for healing.

Patients should be warned that if they decide to explant the product, they should consider the risks of undergoing a second surgical procedure. Implant extraction must be accompanied by an appropriate surgical and post-surgical approach.

Performing additional surgeries does not guarantee a complete repair of the complications.

11. SURGICAL TECHNIQUE

The surgeon must select the appropriate approach based on the repair to be performed. The surgeon may opt for an open or arthroscopic approach, as the surgeon prefers.

Additionally, the surgeon must select the appropriate implant type and dimensions based on the type of injury, the anatomical region, the medical history of the patient and other criteria that the professional may deem important for a satisfactory use of the product during the surgical procedure.

The following surgical technique instructions refer only to product implantation, not considering other procedures that are critical to perform the surgical procedure, such as patient preparation, working portals, procedure finalization, etc.

Some surgical procedures associated to tenodesis PEEK anchors are suggested below. However, other procedures may be performed if the physician deems it appropriate as long as the indications for use are followed:

Shoulder	<ul style="list-style-type: none"> • Proximal biceps tenodesis
Knee	<ul style="list-style-type: none"> • Lateral extra-articular tenodesis (LET) • Posterolateral complex reconstruction • Posteromedial complex reconstruction • Medial patellofemoral ligament (MPFL) reconstruction
Ankle	<ul style="list-style-type: none"> • Tenodesis for flexor <i>hallucis longus</i> (FHL) tendon transfer

Follow these steps for proper product implantation:

For SpinFit PK Tenodesis

- 1) Select the implant insertion site, identify the tendon and prepare it.

NOTE: Using a length of no less than 20 mm of CordFit Infinity #2 (USP standard) (FG-SU040102) suture is recommended to suture the graft, leaving a section of suture available at the end.

- 2) Measure the diameter of the graft on the sutured end using the graft sizing block (FG-IN520020).
- 3) Select the appropriate screw for the graft diameter.
- 4) Using the guide pin (FG-CN305401/FG-CN305402), prepare a hole in the implantation site.
- 5) The hole and the graft position will vary depending on the product model:
 - a. For SpinFit PK Tenodesis with suture eyelet
 - a.i. Without removing the guide pin, select a drill bit with the same diameter as the screw and drill a hole in the bone with a depth of, at least, 23 mm.
 - a.ii. Thread the graft preparation sutures through the eyelet of the device so that the eyelet is positioned at the end of the graft, holding the graft preparation sutures. Tension the graft sutures to prevent the graft from moving.
 - a.iii. Uncouple the eyelet retention suture from the proximal end of the handle and traction until the eyelet is closed, in order to secure the end of the graft tightly.

a.iv. Position the insertion device together with the end of the graft inside the hole in the bone until the laser marking is reached while making sure that the eyelet suture tension is maintained.

ATTENTION: Adjust graft tension at this point, since once the implant is placed, the tension cannot be modified.

b. For SpinFit PK Tenodesis with PEEK eyelet

b.i. Without removing the guide pin, select a drill bit with the same diameter as the screw and drill a hole in the bone with a depth of, at least, 25 mm.

b.ii. Place the graft using the fork eyelet to insert it up to the laser marking. **NOTE:** Graft preparation sutures may be threaded in the side holes of the PEEK eyelet.

ATTENTION: Adjust graft tension at this point, since once the implant is placed, the tension cannot be modified.

b.iii. Uncouple the eyelet retention suture from the proximal end of the handle.

- 6) Turn the proximal end of the insertion device clockwise while holding the handle to allow for the implant to move into the hole in the bone until it has been completely inserted.
- 7) Take the insertion device out and remove the excess suture.

After use, discard the insertion device, excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

For Spirom PK Tenodesis

- 1) Select the implant insertion site, identify the tendon and prepare it.
NOTE: Using a length of no less than 20 mm of CordFit Infinity #2 (USP standard) (FG-SU040102) suture is recommended to suture the graft, leaving a section of suture available at the end.
- 2) Measure the diameter of the graft on the sutured end using the graft sizing block (FG-IN520020).
- 3) Select the appropriate screw for the graft diameter.
- 4) Using the guide pin (FG-CN305401/FG-CN305402), prepare a hole in the implantation site.
- 5) Without removing the guide pin, select a drill bit with the same diameter as the screw and drill a hole in the bone with a depth of, at least, 23 mm.

- 6) Attach the implant to the Fergus hex screwdriver (FG-IN601206) and position the graft inside the hole in the bone.

ATTENTION: Adjust graft tension at this point, since once the implant is placed, the tension cannot be modified.

- 7) Turn the instrument clockwise to advance the implant into the hole in the bone until it has been completely inserted.
- 8) Take the instrument out and remove the excess suture.

After the procedure, discard the excess sutures, contaminated disposable surgical material and packaging per local hospital, administrative or governmental procedures.

12. MATERIAL SPECIFICATIONS

All sizes of Spirom PK Tenodesis and SpinFit PK Tenodesis screws are made with polyetheretherketone (PEEK) as per the ASTM F2026 standard.

The eyelet of SpinFit PK Tenodesis may be made of PEEK (fork eyelet) or CordFit #2 suture (USP standard), which is manufactured using ultra-high-molecular-weight polyethylene (UHMWPE)

The UHMWPE suture comes in white-blue. The *Cobalt Chromite Blue-Green Spinel* pigment is used for the blue threads. White threads are not pigmented.

The insertion device is made up of the following materials:

- Handle: it is manufactured with medical grade ABS. There are different colors to identify the implant diameter for which the insertion device is used. The colors used to identify the products are mentioned below:
 - The 6 mm SpinFit PK Tenodesis model comes in gray.
 - The 7 mm SpinFit PK Tenodesis model comes in light blue.
 - The 8 mm SpinFit PK Tenodesis model comes in blue.
- Capillary: manufactured with stainless steel, as per the AISI 304 standard.

Check the packaging labels for the implant materials.

13. LEARN MORE

To learn more about this product, visit our website: www.fergusmedical.com

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